



3-layer Protective Mask

Executive standard:GB2626-2006

Comply standard:YY0469-2011

External side



Internal side



Size:175mm(L)×95mm(W)

Packing:50pcs/box

Material: Layer 1&3: Polypropylene spunbond Layer 2: Meltblown cloth

Executive standard:GB2626-2006 Comply standard:YY0469-2011

Testing Result :

*Under testing report No. SHF19080222-02R1, by ICAS TESTING TECHNOLOGY SERVICE (SHANGHAI) CO., LTD.

No.	Item	Standard Clause	Technical Requirement	Testing Result	Conclusion
1	Inhalation resistance	2.2.3	≤175	32.4	Pass
2	Expiratory resistance	2.2.3	≤145	16.3	Pass
3	Strength of mask band and connection of mask band and mask body	2.2.3	≥20	31	Pass
4	View under mask	2.2.3	≥60°	65°	Pass
5	Filtration efficiency	2.2.4	Lithologic medium Oily medium	99.9 99.9	Pass
6	Protective effect	2.2.5	≥90%	95	Pass
7	PH	2.2.6	4.0-4.8	4.5	Pass
8	Label	3	Manufacturer name, date of production, QC certificate	Comply	Pass

*Gamma Ray Sterilized

3-layer construction



Leakage test by water



Vacuum package in bag



Packed in box



Production is finished in cleanness grade workshop

