

LabSafe Series

-25°C/Under-counter/ MDF-25U150

ANTECH undercounter lab freezer MDF-25U150 is ideal for small medical practices, pharmacies, and labs. Its height is designed to match laboratory table and optimized usage of limited space. Despite a compact freezer, the functions are complete for medical and laboratory applications, including data storage and USB port for download, high/low temperature alarm, door open alarm, power failure alarm, remote alarm contact, etc. Therefore, MDF-25U150 meets the most rigorous requirements of healthcare and scientific research.

This unit uses environmentally friendly, natural refrigerants that are fully compliant with new EPA/SNAP/UL regulations.

MDF-25U150 is a stackable freezer (using Antech stacking kits) to install in a limited space. You can order with various combinations with flexibility for your specific demand.

▼ MDF-25U150



▼ With Stacking Kits



NATURAL
REFRIGERANT



USB PORT



HIGH/LOW
TEMPERATURE ALARM



DOOR AJAR A
LARM



REVERSIBLE
DOORS



STACKABLE
DESIGN

MODEL	MDF-25U150
Temperature range (°C)	-10~-25
Capacity (L)	120
External dimensions W*D*H (mm)	600×560×820
Shelves (pcs)	3
Data storage & USB port	Yes
Power supply	220V/50Hz; 220V/60Hz; 120V/60Hz