

# LabSafe Series

## -25°C/ Desktop/ MDF-25U60

ANTECH desktop lab freezer MDF-25U60 is ideal for small medical practices, pharmacies, and labs. It is a great choice when you need a small lab freezer for specialized medicines, vaccines or lab materials. It is perfectly suited for additional storage where space is a consideration. The MDF-25U60 is designed for the most rigorous life science, research, healthcare and clinical applications.

This unit uses environmentally friendly, natural refrigerants that are fully compliant with new EPA/SNAP/UL regulations.



▶ MDF-25U60



▶ With Stacking Kits



CFC /HCFC  
FREE



USB PORT



POWER FAILURE  
ALARM



DOOR AJAR A  
LARM



REVERSIBLE  
DOORS



STACKABLE  
DESIGN

MODEL	MDF-25U60
Temperature range (°C)	-10~-25
Capacity (L)	65
External dimensions W*D*H (mm)	510×560×670
Shelves (pcs)	2
Data storage & USB port	Yes
Power supply	220V/50Hz; 220V/60Hz; 120V/60Hz