


Racks for Antech Freezers

Model	Rack type (boxes each layer x layers)	Rack model		Racks per freezer	2" box storage capacity	Maximum 2ml vial capacity
		Sliding drawer	Side access			
MDF-86U108	3x5	VFD-352	VFF-352	4	60	6,000
MDF-86U158	3x4	VFD-342	VFF-342	9	108	10,800
MDF-86U340T	4x4	VFD-442	VFF-442	12	192	19,200
MDF-86U458T/458D	4x5 & 4x4	VFD-452&VFD-442	VFF-452&VFF-442	8+8	288	28,800
MDF-86U588T/588D	5x5	VFD-552	VFF-552	16	400	40,000
MDF-86U728T/728D	4x5 & 4x6	VFD-452&VFD-462	VFF-452&VFF-462	18+6	504	50,400
MDF-86U838T/838D	5x5	VFD-552	VFF-552	24	600	60,000
MDF-86C118	1x8	N/A	CFC-8-2	12	96	9,600
MDF-86C458	1x10	N/A	CFC-10-2	36	360	36,000
MDF-86C728	1x13	N/A	CFC-13-2	45	585	58,500
MDF-150C118	1x9	N/A	CFC-9-2	6	54	5,400
MDF-150C208	1x11	N/A	CFC-11-2	15	165	16,500

Vertical Racks with Spring Clip for 2" Boxes



- ▶ Spring clips for easy operation and prevent box drop.
- ▶ Dual large handles at top for easy handling.
- ▶ Label holders for easy identification.
- ▶ Constructed of high-grade corrosion resistant stainless steel.
- ▶ Precision craftsmanship & no sharp edges.

Vertical Racks for 2" Boxes

Item#	Description / Capacity	H x W x D, mm
CFC-7-2	7 Boxes	432 x 140 x 140
CFC-8-2	8 Boxes	488 x 140 x 140
CFC-9-2	9 Boxes	544 x 140 x 140
CFC-10-2	10 Boxes	600 x 140 x 140
CFC-11-2	11 Boxes	656 x 140 x 140
CFC-12-2	12 Boxes	712 x 140 x 140
CFC-13-2	13 Boxes	768 x 140 x 140
CFC-14-2	14 Boxes	824 x 140 x 140
CFC-15-2	15 Boxes	880 x 140 x 140
CFC-16-2	16 Boxes	936 x 140 x 140

Note: including handle height 40mm.

Sliding Drawer Freezer Racks for 2" & 3" Boxes



- ▶ Handle at each drawer for easy handling.
- ▶ Label holder at each drawer for easy identification.
- ▶ Constructed of high-grade corrosion resistant stainless steel.
- ▶ Precision craftsmanship & no sharp edges.

For 2" Freezer Boxes

Item#	Boxes per layer x Layers	L x H x W, mm
VFD-332	3 x 3 Rack/ for 9 boxes	424 x 177 x 140
VFD-342	3 x 4 Rack/ for 12 boxes	424 x 236 x 140
VFD-352	3 x 5 Rack/ for 15 boxes	424 x 295 x 140
VFD-432	4 x 3 Rack/ for 12 boxes	557 x 177 x 140
VFD-442	4 x 4 Rack/ for 16 boxes	557 x 236 x 140
VFD-452	4 x 5 Rack/ for 20 boxes	557 x 295 x 140
VFD-462	4 x 6 Rack/ for 24 boxes	557 x 354 x 140
VFD-542	5 x 4 Rack/ for 20 boxes	690 x 236 x 140
VFD-552	5 x 5 Rack/ for 25 boxes	690 x 295 x 140
VFD-652	6 x 5 Rack/ for 30 boxes	823 x 295 x 140

For 3.75" Freezer Boxes

Item#	Boxes per layer x Layers	L x H x W, mm
VFD-323	3 x 2 Rack/ for 6 boxes	424 x 202 x 140
VFD-333	3 x 3 Rack/ for 9 boxes	424 x 303 x 140
VFD-423	4 x 2 Rack/ for 8 boxes	557 x 202 x 140
VFD-433	4 x 3 Rack/ for 12 boxes	557 x 303 x 140
VFD-443	4 x 4 Rack/ for 16 boxes	557 x 404 x 140
VFD-453	4 x 5 Rack/ for 20 boxes	557 x 505 x 140
VFD-523	5 x 2 Rack/ for 10 boxes	690 x 202 x 140
VFD-533	5 x 3 Rack/ for 15 boxes	690 x 303 x 140

Side Access Freezer (Frame Type) Racks for 2" & 3" Boxes



- ▶ Large handles for easy handling.
- ▶ Label holders for easy identification.
- ▶ Constructed of high-grade corrosion resistant stainless steel.
- ▶ Precision craftsmanship & no sharp edges.

For 2" Freezer Boxes

Item#	Boxes per layer x Layers	L x H x W, mm
VFF-332	3 x 3 Rack/ for 9 boxes	425 x 171 x 137
VFF-342	3 x 4 Rack/ for 12 boxes	422 x 229 x 137
VFF-352	3 x 5 Rack/ for 15 boxes	422 x 284 x 137
VFF-432	4 x 3 Rack/ for 12 boxes	560 x 171 x 137
VFF-442	4 x 4 Rack/ for 16 boxes	560 x 229 x 137
VFF-452	4 x 5 Rack/ for 20 boxes	560 x 284 x 137
VFF-462	4 x 6 Rack/ for 24 boxes	560 x 339 x 137
VFF-542	5 x 4 Rack/ for 20 boxes	696 x 229 x 137
VFF-552	5 x 5 Rack/ for 25 boxes	696 x 284 x 137
VFF-652	5 x 6 Rack/ for 30 boxes	696 x 339 x 137

For 3.75" Freezer Boxes

Item#	Boxes per layer x Layers	L x H x W, mm
VFF-323	3 x 2 Rack/ for 6 boxes	424 x 202 x 140
VFF-333	3 x 3 Rack/ for 9 boxes	424 x 303 x 140
VFF-423	4 x 2 Rack/ for 8 boxes	557 x 202 x 140
VFF-433	4 x 3 Rack/ for 12 boxes	557 x 303 x 140
VFF-443	4 x 4 Rack/ for 16 boxes	557 x 404 x 140
VFF-453	4 x 5 Rack/ for 20 boxes	557 x 505 x 140
VFF-523	5 x 2 Rack/ for 10 boxes	690 x 202 x 140
VFF-533	5 x 3 Rack/ for 15 boxes	690 x 303 x 140